


Top Doc Magnetic Surgery Game

Carefully remove game table from the delivery box and extend & lock the table legs.



Attached the velcro skirt around the table. Plug in the game to the outlet and attached the wand to the game (towards the head of the patient on the side of the table).



You are ready to play!



Easy Pieces are engraved with E

Want more of a challenge - Try the Medium Pieces engraved with M



Top Doc takes a classic board game, and scales it up while adding a technology twist for greater appeal and improved game play. With Top Doc, the players wand is actually an electromagnet. So when the player touches the wand to a game piece, it automatically latches onto the wand and magnetizes. While carefully removing the piece, should contact be made with the cutout, the magnet automatically turns off and the piece drops. An audible beep is triggered along with a visual light, and the unit records the mistake on a large 3 inch LED display. Games are available in Easy and Medium difficulty levels, and can be easily changed by swapping out the game pieces.

Contents:

- 1 - Carry Case
- 1 - 4 foot Game Table
- 1 - Table Skirt
- 1 - Set of 9 Easy Pieces
- 1 - Set of 9 Medium Pieces
- 1 - Red Zipper Pouch for 18 Game Pieces
- 1 - Power Cord
- 1 - Magnetic Wand
- 1 - Blue Zipper Pouch for Cord & Wand

